

HOW TO USE

AT-HOME COVID-19 RAPID TESTS

At-home rapid tests can help inform our decisions, so we can avoid spreading COVID-19, and can get care ASAP.

Rapid tests, aka antigen tests, produce results after 15 minutes, and are approved for ages 2+.

If possible, take a second test in a few days, because rapid tests aren't quite as accurate as PCR tests.

1

Make sure that the seal is not broken. Open the package up, and read the directions.

2

Lay out all the equipment, labeled to the right. First, open the cartridge and lay it flat.

3

Use the applicator to put the designated amount of drops on the marked cartridge section.

4

Open the swab's container and insert the swab about $\frac{3}{4}$ " up each nostril. In each nostril, rotate the swab a few times, then insert it into the cartridge.

5

Seal your test, following the included directions. After time is up, check the cartridge and reference the included key to determine your results.

cartridge

swab

applicator

